 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
ProView[®] Minimal Access Portal (MAP) System
Tubular Retractor System
Expandable Retractor System
Percutaneous Screw Delivery System
ONYX[™] Instrumentation Set

 Click directory below for desired language

English 	2-3
Español 	4-5
Italiano 	6-7
Deutsch 	8-9
Português 	10-11
日本語 	12-13
Tiếng Việt 	14-15

INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

ProView® Minimal Access Portal (MAP) System Tubular Retractor System Expandable Retractor System Percutaneous Screw Delivery System ONYX™ Instrumentation Set

Description:

The ProView MAP Systems consist of surgical instruments intended to aid surgeons in visualization of surgical sites for spinal procedures.

The ProView Tubular Retractor System includes stainless steel tubular retractors in multiple lengths and diameters, a handle for insertion into the operative site, fiber optic lighting, and table attachments to connect to the side rail of the operating room table.

The ProView Expandable Retractor System includes detachable blades in multiple lengths, fiber optic lighting, and table attachments to connect to the side rail of the operating room table.

The ProView Percutaneous Screw Delivery System consists of instrumentation designed to deliver multi-axial pedicle screws and rods in a minimally invasive manner. The system is designed to reduce muscle trauma compared to the traditional open surgical procedure.

The ProView ONYX System features a comprehensive set of black-coated Kerrisons, curettes, pituitaries, stylets, hand held retractors, probes and scalpel handles. The instruments are designed to perform specific functions such as cutting, grasping, dissecting, probing, retracting and draining.

Indications for Use:

The ProView Tubular and Expandable Retractor System are intended to aid the surgeon's visualization of the surgical area and allow for performance of spinal procedures such as herniated disc repair, visualization of the circumferential decompression of the nerve roots, aiding in the search and removal of nucleus material, spinal fusion, or insertion of spinal implants.

The ProView Percutaneous Screw Delivery System is used to deliver multi-axial pedicle screws and rods in a minimally invasive manner.

The ProView ONYX System is intended to aid in cutting and removal of soft tissue during a surgical procedure.

Contraindications:

The ProView MAP System is contraindicated for use in patients with:

1. Morbid obesity.
2. Mental illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia.
7. Patients unwilling or unable to follow post-operative care instructions.
8. Any circumstances not listed under the heading Indications for Use.

Potential Adverse Events:

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Neurological injury.
2. Vascular or visceral injury.
3. Foreign body (allergic) reaction to instruments, debris, etc.
4. Infection.
5. Hemorrhage.
6. Death.

Note: Potential risks identified with the use of the device system may require additional surgery.

Warnings and Precautions:

1. The ProView MAP System is sold nonsterile and therefore must be thoroughly cleaned and sterilized before each use.
2. Care should be exercised in the handling and storage of instruments. Instruments should not be scratched, notched, or otherwise damaged since such actions may reduce functional performance. Store away from corrosive environments.
3. The retractors should be assembled prior to surgery.
4. If used around the spinal cord and nerve roots, extreme caution should be taken.
5. The reusable fiber optic light cables are designed for use with 300 watt xenon illuminators using the provided ACMI adapter cable. Do not use light sources rated higher than 300 watts or any cables other than the provided ACMI adapter cable. Use of higher watt sources or cables other than the provided ACMI adapter cable could result in overheating, causing product failure and patient injury.
6. Do not operate the light source and adapter cable without the reusable fiber optic light cables attached. Without the reusable fiber optic light cable, the output from the adapter cable is extremely bright, hot, and may cause burns, ignite drapes/gowns, or temporarily blind vision.

Cleaning:

All instruments must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse.

Whenever possible, do not allow the blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use. The use of an instrument lubricant, that is compatible with the method of sterilization to be used, is recommended before instruments are sterilized. Proper application of lubricants is an integral step in maintaining the long life of the instrument.

Instruments should never be exposed to cleaning agents containing any peroxides. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

Sterilization:

The ProView MAP System is supplied NON-STERILE. Prior to use, all instruments should be placed in the instrument case which will be wrapped in an FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The ProView MAP System instruments are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.

Latex Information:

The instruments and/or packaging material for the ProView MAP System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

Rx Only		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Manufacturer
	Orthofix.com/IFU		Authorized Representative
	Catalogue Number		Lot Number
	Provided Non-Sterile		

 Orthofix Inc.
 3451 Plano Parkway
 Lewisville, Texas 75056-9453 U.S.A.
 1-214-937-3199
 1-888-298-5700
 www.orthofix.com

Australian Sponsor
 Emergo Australia
 Level 20, Tower II
 Darling Park
 201 Sussex Street
 Sydney, NSW 2000
 Australia

 Medical Device Safety Service
 (MDSS)
 Schiffgraben 41
 30175 Hannover
 Germany
 +49 511 6262 8630
 www.mdss.com

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