Orthofix
3451 Plano Parkway
Lewisville, Texas 75056-9453 USA
214-937-2000
1-888-298-5700
www.orthofix.com

EC REP Medical Device Safety Service (MDSS)
Schiffgraben 41
30175, Hannover
Germany

Rx Only
CE 0086

English

Sterilization in Rigid Sterilization Containers

Orthofix implants and instruments are supplied NON-STERILE and must be thoroughly cleaned prior to each use. Prior to sterilization, all implants and instruments should be placed in their instrumentation / implant case.

Rigid sterilization container must be cleaned, inspected and prepared according to the rigid sterilization manufacturer instructions for use.

Orthofix system is validated for sterilization in rigid sterilization containers. Select the appropriate rigid sterilization container (either Filtered Bottom or Solid Bottom) to properly enclose the Orthofix instrumentation case(s) (recommended 23 1/4" long x 11 1/4" wide container). The following sterilization cycle was validated:

Sterilization Method — Steam Prevac Cycle
Temperature: 270° F (132° C)
Preconditioning: 4 Pulses
Exposure time: 4 minutes
Dry Time: 30 minutes

Note: Rigid Sterilization Containers with Solid Bottom cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed as per ANSI/AAMI ST79: "Comprehensive guide to steam sterilization and sterility assurance in health care facilities". Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Note: Please refer to the "Instructions for Use" supplied with the Orthofix system for full information on indications for use, contraindications, warnings, precautions, adverse reactions information, sterilization and other useful information.

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