



# IsoTis<sup>®</sup> Pure Strip

Demineralized Bone Matrix Strip



IsoTis<sup>®</sup> Pure Strip

# IsoTis<sup>®</sup> Pure Strip

## DESIGN RATIONALE

IsoTis<sup>®</sup> Pure Strips are scaffolds with verified osteoinductive\* potential. Composed of 100% human tissue, the implants are designed specifically for use as bone void fillers in lumbar fusion procedures.

## Demineralized Bone Matrix Strip

## FEATURES

### Flexible Handling

- The hydrated graft is pliable and may be molded or cut as desired
- Conforms to the defect site and make contact with the decorticated bone
- The graft is cohesive and maintains integrity during irrigation


### Hydration

- Hydrate with sterile saline, sterile water, blood, bone marrow aspirate (BMA) or platelet rich plasma (PRP)
- Hydrating fluid is retained within the implant at the graft site
- Graft hydrates in at least 3 minutes
- Alternative hydration options provide flexibility
  - May be hydrated *in situ*
  - Shorter hydration times will result in firmer handling

### Verified Osteoinductive Potential

- This 100% human tissue is a scaffold with verified osteoinductive\* potential
- The Demineralized Bone Matrix (DBM) is sterilized through the Cancell<sup>®</sup> SP DBM sterilization process, which is designed to preserve naturally occurring protein and physical structure of the tissue
- Osteoinductive (OI) potential is verified by 100% lot testing

\*DBM or representative finished implant is either assayed *in vivo* in the modified athymic nude rat for bone formation or *in vitro* for endogenous BMP-2 as a surrogate test marker for osteoinductive potential. Because the combination of various proteins is responsible for osteoinductive potential, DBM when assayed *in vitro*, is also screened for the presence of BMP-7. Findings from an *in vitro* assay or animal model are not necessarily predictive of human clinical results.




## SAFETY

The highest level of safety is provided through redundant safeguards, achieving a Sterility Assurance Level (SAL) of 10<sup>-6</sup>:

- Stringent donor screening
- Independent CLIA Laboratory Testing
- BioCleanse<sup>®</sup> Tissue Sterilization Process
- Cancell<sup>®</sup> SP DBM Sterilization Process

## Sizing

- Use 6cc size for filling bone voids in single level fusion procedures
- Use the 12cc or 8cc x 2 unit configurations for filling bone voids for longer fusion or deformity procedures



Due to the inherent variation of human tissue, actual dimensions may vary.

### IsoTis Pure Strip

Reference	Description
06-2000-060	IsoTis® Pure Strip, 6cc (approx 5cm x 2.5cm)
06-2000-120	IsoTis® Pure Strip, 12cc (approx 10cm x 2.5cm)
06-2000-160	IsoTis® Pure Strip, 16cc (approx 10cm x 1 cm) x 2 units



For more information or to place an order, please contact:

IsoTis OrthoBiologics, Inc 2 Goodyear, Suite A, Irvine CA 92618

Phone 800.550.7155

[SeaSpine.com](http://SeaSpine.com)

BioCleanse and Cancellé are registered trademarks of RTI Surgical, Inc. IsoTis is a registered trademark of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. The SeaSpine logo and the IsoTis logo are trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2016 SeaSpine Orthopedics Corporation. All rights reserved. D0000445A 2016-02