INSTRUCTIONS FOR USE
Important Information - Please Read Prior to Use

Rx Only
CE 0086

 Orthofix US LLC
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
**CONSTRUX® Mini PEEK Vertebral Body
Replacement System**

 Click directory below for desired language

English	EN	2-3
Español	ES	4-5
Italiano	IT	6-7
Deutsch	DE	8-9
Português	PT	10-11
英文	ZH-TW	12-13

INSTRUCTIONS FOR USE

Important Information - Please Read Prior to Use

Rx Only
CE 0086

Orthofix US LLC
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

EC REP Medical Device Safety Service (MDSS)
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

CONSTRUX® Mini PEEK Vertebral Body Replacement System

Description:

The CONSTRUX Mini PEEK VBR (Vertebral Body Replacement) System is comprised of a variety of implants manufactured from PEEK (Polyetheretherketone), as described by ASTM F2026, with titanium markers as described by ASTM F67. The implants are available in two footprint sizes, a small and a large. The implants are available in various heights, in one-millimeter increments. The superior and inferior surfaces of the implant have a pattern of ripples to provide increased stability and help prevent anterior/posterior movement of the device.

The CONSTRUX Mini PEEK VBR System is intended for vertebral body replacement to aid in the surgical correction and stabilization of the spine.

The CONSTRUX Mini PEEK VBR System is not intended to be used as a stand-alone device. The system must be used with supplemental fixation. The CONSTRUX Mini PEEK VBR device is used singly and is implanted using an anterior approach. The system is provided non-sterile and requires sterilization prior to use.

Indications for Use:

The CONSTRUX Mini PEEK VBR System is indicated for use in the thoracolumbar spine (T1-L5) for partial replacement (i.e., partial vertebrectomy) of a diseased vertebral body resected or excised for the treatment of tumors, to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body. The system is also indicated for treating fractures of the thoracic and lumbar spine.

The CONSTRUX Mini PEEK VBR System is designed to restore the biomechanical integrity of the anterior, middle, and posterior spinal column even in the absence of fusion for a prolonged period of time. The Partial VBR Device is intended to be used with autograft or allograft.

The CONSTRUX Mini PEEK VBR System is intended for use with supplemental fixation. The supplemental fixation system that may be used with the CONSTRUX Mini PEEK VBR System is the Orthofix Spinal Fixation System (SFS) or the Firebird Spinal Fixation System.

Contraindications:

The CONSTRUX Mini PEEK VBR System, as with other orthopedic implants, is contraindicated for use in patients with:

1. Active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection.
2. Rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis which may prevent adequate fixation.
3. Conditions that may place excessive stresses on bone and implants, such as severe obesity, pregnancy or degenerative diseases. The decision to use this system in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
4. Prior fusion at the level to be treated.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

Note: As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, permanent pain and/or deformity. Rarely, some complications may be fatal.

Note: The CONSTRUX Mini PEEK VBR System is intended for use in affected vertebral body segments that are equal to or smaller than the size of the device. For larger affected vertebral body segments, a larger device indicated for partial or full VBR is recommended.

Warnings and Precautions:

The surgeon should be aware of the following when using implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. The size and shape of human bones present limiting restrictions on the size and strength of implants. No implant can be expected to withstand the unsupported stresses of full weight bearing.
2. Correct handling of the implant is extremely important. Implants should not be bent, notched or scratched. These operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.
3. Single use only. No surgical implants should be reused. Any implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns that may lead to fatigue failure.
4. Non-sterile; the CONSTRUX Mini PEEK VBR System implants and instruments are provided non-sterile, and therefore, must be sterilized before each use.
5. Postoperative care is important. The patient should be instructed in the limitations of the implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
6. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.

MRI Compatibility Information:

The CONSTRUX Mini PEEK VBR System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

Cleaning:

CONSTRUX Mini PEEK VBR System implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluids, it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned and sterilized after each use. Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning process described below.

None of the instruments require disassembly prior to cleaning.

From Point of Use:

Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Preparation for Cleaning:

1. Any instruments with moving parts (e.g., knobs, triggers and hinges) should be separated or placed in the open position to allow better access of the cleaning fluid to the difficult to clean areas.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or plastic bristle brush to remove any visible soil from the instruments.

Manual Cleaning:

1. Because none of the instruments can be disassembled, use a soft cloth and/or a soft plastic bristle brush to remove any visible soil from the outside and inside of the instrument.
2. Prepare Vesphene® Ise agent at the use-dilution recommended by label direction – 1 ounce per gallon or 1 mL of Vesphene® Ise to 128 mL of potable tap water per manufacturer's recommendations.
3. Bathe instruments in prepared room temperature solution as recommended by the detergent manufacturer.
4. Manually agitate instruments in Vesphene® solution for 1.5 minutes.
5. Scrub instruments with a soft plastic bristle brush if visible soil is noted and use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the inferior lumen to remove the soil.
6. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
7. Hang dry the device.
8. Visually inspect the instruments for visible soil.
9. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. Use a soft cloth or soft plastic bristle brush to remove any visible soil from the outside and inside of the instruments. Use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil from the inner lumens prior to automated cleaning.
2. Set up the washer / disinfectant detergent dose as ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer recommendations.
3. Remove the instruments from the instrument case and place them into the washer baskets.
4. Orient instruments in the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
 - a. Pre-Wash 1: cold potable water, 2 minutes
 - b. Enzyme/Detergent treatment:
 1. Spray, 20 seconds
 2. Soak, 1 minute
 3. Rinse cold potable water, 15 seconds
 4. Rinse cold potable water, 15 seconds
 - c. Wash ≥ 65°C, 2 minutes using Endozime AW Plus®
 - d. Rinse 1: hot potable water, 15 seconds
 - e. Rinse 2: hot potable water, 15 seconds
 - f. Rinse 3: hot potable water, 15 seconds
 - g. Rinse 4: hot potable water, 15 seconds
 - h. Thermal rinse ≥ 93°C, 1 minute
 - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
 - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
 - k. Dry at 115°C, 7 minutes
6. Visually inspect the instruments for visible soil
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix if instruments are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

Sterilization:

Sterilization in Orthofix Cases with Blue Wrap:

The CONSTRUX Mini PEEK VBR System implants and instruments are supplied NON-STERILE. Prior to use, all implants and instruments should be placed in the instrument/implant case which will be wrapped in an FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

- | | | |
|----------------------------|-----|----------------------------|
| Method: Steam | or: | Method: Steam |
| Cycle: Gravity | | Cycle: Prevac |
| Temperature: 270°F (132°C) | | Temperature: 270°F (132°C) |
| Exposure time: 15 minutes | | Preconditioning: 4 Pulses |
| Drying time: 30 minutes | | Exposure time: 4 minutes |
| Double wrapped | | Drying time: 30 minutes |
| | | Double wrapped |

Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer's instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix instrumentation cases (recommended 23¼" long x 11¼" wide container). The following sterilization cycle has been validated:

- Method: Steam
 Cycle: Prevac
 Temperature: 270°F (132°C)
 Preconditioning: 4 Pulses
 Exposure time: 4 minutes
 Drying time: 30 minutes

Note: Rigid Sterilization Containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed as per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The CONSTRUX Mini PEEK VBR System instruments and implants are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Healthcare Professional (e.g., customer or user of this system of products) who has any complaints, or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix US LLC, 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by email at complaints@orthofix.com.

Further Information:

A recommended operative technique for the use of this system is available upon request. Please contact Orthofix using the information listed above.

Latex Information:

The implants, instruments and/or packaging material for the CONSTRUX Mini PEEK VBR System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Rx Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician		
	See Instructions for Use	REF	Catalogue Number
	Orthofix.com/IFU		Manufacturer
	Single Use Only Do Not Reuse	EC REP	Authorized Representative
	Provided Non-Sterile	SN	Serial Number
LOT	Lot Number		

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