INSTRUCTIONS FOR USE

Important Information — Please Read Prior to Use

Orthofix Inc.

3451 Plains Parkway

Lewisville, Texas 75056-9453 USA

214-937-2000

complaints@orthofix.com

Medical Device Safety Service (MDSS)

Schilligstrasse 11

30175, Hannover

Germany

Rx Only

0086

English

EN

Device System Name:
BlackJack™ Anterior Lumbar Distractor/Inserter System

Description: The BlackJack instrument integrates methods of distraction and insertion into one easy-to-use tool for placement of Intervertebral Body Devices, (IBD) partial Vertebral Body Replacement (VBR) implants or allograft products during anterior lumbar procedures. By simply turning the handle of the BlackJack instrument, the product is advanced down the track, acting as a wedge, distracting the blades and the selected spinal space as dictated by the height and lordotic angle of the product. It then gently places the product into position and removes the blades from the selected space allowing the bone to retract back onto the implanted product, securing it in place.

The BlackJack instrument system includes one distraction/insertion instrument, consisting of a handle, stainless steel blades and threading rod with a unique locking mechanism to prevent the blades from opening during transfer to the surgeon. Each system also contains one tip for AlloQuent allograft products that will fit all sizes and three tips to fit the various footprints for PILLAR™ AL and SA implants.

Refer to package Instructions for Use and inserts provided with the specific implantable allograft, IBD,partial VBR products for information specific to the product, including: indications for use, contraindications, warnings, precautions, adverse reaction information and sterilization information.

Contraindications:
Contraindications of the BlackJack Distractor/Inserter System include:

1. Instruments should not be used for anything other than their intended use
2. Morbid obesity
3. Mental illness
4. Alcoholism or drug abuse
5. Pregnancy
6. Metal sensitivity/allergies
7. Severe osteopenia
8. Patients unwilling or unable to follow post-operative care instructions
9. Active infection
10. Any circumstances not listed under the heading Indications.

Potential Adverse Effects:
As with any major surgical procedure involving spinal instrumentation, infrequent possible adverse effects include, but are not limited to:

1. Device component fracture
2. Neurological injury
3. Vascular or visceral injury
4. Foreign body (allergic) reaction to instruments, debris, corrosion products, including metallosis, straining, tumor formation, and/or autoimmune disease
5. Infection can result if instruments are not properly cleaned and sterilized
6. Hemorrhage
7. Cessation of any potential growth of the operated portion of the spine
8. Death
9. Potential complications associated with the implanted product are described in the Package Insert or Instructions for Use provided with the implantable product.

Warnings and Precautions:
The BlackJack Distractor/Inserter System should only be used by surgeons who are experienced in the surgical procedure and have undergone adequate training with the device. The surgeon should be aware of the following when using the instrument:

1. The BlackJack Anterior Lumbar Distractor/Inserter System instruments must be cleaned and sterilized prior to use
2. Use the instruments as intended and not under impaction unless specifically designed for this purpose
3. If used around the spinal cord and nerve roots, extreme caution should be taken
4. Care should be exercised in the handling and storage of instruments. Instruments should not be scratched, notched, or otherwise damaged since such actions may reduce functional performance. Store away from corrosive environments.
5. Instruments should be carefully placed in the cases and stored in a cool dry place when not in use.
6. Instruments should be examined for functionality or damage prior to use and any damaged instruments should not be used to implant devices.
7. Improper handling of instruments may result in malfunction during surgery.
8. Damage to implants from instruments must be avoided, or device failure could occur.
9. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not attempt to re-sterilize single-use implants that come in contact with body fluids.

Cleaning:
All instruments and implants must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments and implants that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

Sterilization: The BlackJack Anterior Lumbar Distractor/Inserter System components are supplied NON-STERILE. Prior to use, all components should be steam sterilized by the hospital using the recommended cycle:

Method: Steam	Or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250° F (121° C)		Temperature: 270° F (132° C)
Exposure time: 30 minutes		Exposure time: 8 minutes

Product Complaints: Any Healthcare Professional (e.g., customer or user of this system), who has any complaints, or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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SI-1024.20-PL-US © Orthofix Inc.
RX Only Federal (U.S.A.) law restricts this device to sale by or on the order of a physician