

Fitting & Care Instructions - Flex TLSO SP

Model: DCT-64

Instructions for putting on brace



Step 1: Open brace and lay flat on a table. Detach side wings from right and left to find your waist circumference and size.



Step 2: Place your head through opening and have the front T-bar on your abdomen and chest.



Step 3: Wrap left side wing around waist and Velcro closed on top of the abdominal panel.



Step 4: Locate the rubber pull tabs on either side of the brace and pull to desired support.



Step 5: Open the back flap of the brace and locate the Velcro strip. Unhook the strip and lower or raise the height of the shoulder straps to the desired level.



Step 6: Locate the bottom and top screws and raise or lower the T-bar to the desired height.



Step 7: Locate the thumb screw on the side of the T-bar and adjust to the desired support.



Step 8: Locate the front and back straps on the top of the shoulder straps and tighten to the desired support.



Size 0	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7
22" - 32"	31" - 37"	36" - 40"	39" - 43"	42" - 46"	45" - 49"	48" - 52"	51" - 55"

Care Instructions: Hand wash with a light detergent and let air dry.

Warranty: 6 month warranty from any manufacturer defects

