INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only

Orthofix US LLC
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

Device System Name:

NGage® Surgical Mesh System

Click directory below for desired language

English EN 2-3

Español ES 4-5

Italiano IT 6-7

Deutsch DE 8-9

Português PT 10-11

Français FR 12-13

AW-74-9902 Rev.AF
NG-1901-PL-US © Orthofix US LLC 10/2026

INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only

Orthofix US LLC
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

English EN

Device System Name:

NGage® Surgical Mesh System

Description:
The NGage Surgical Mesh System is a diamond pattern, surgical mesh device. The body is manufactured from commercially pure (CP) titanium conforming to ASTM F67, the end rings, the standard ring, and the screws are manufactured from titanium alloy that conforms to ASTM F136. Because of the construction, the angle and the length of the mesh can be reduced incrementally to adjust it to individual anatomical conditions. The system is sold non-sterile.

Indications for Use:
The NGage Surgical Mesh System is indicated for use in the thoraco-lumbar spine (T1-L5) to replace a diseased vertebral body resected or excised for the treatment of tumors, to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body.

The NGage Surgical Mesh System is also indicated for treating fractures of the thoracic and lumbar spine.

The NGage Surgical Mesh System is designed to restore the biomechanical integrity of the anterior, middle, and posterior spinal column even in the absence of fusion for a prolonged period. It is recommended to pack bone graft material inside the mesh cage prior to implantation.

The NGage Surgical Mesh System is intended for use with supplemental fixation. As an example, the supplemental fixation system that may be used with the system is the Orthofix Spinal Fixation System.

Contraindications:
The NGage Surgical Mesh System, as with other metallic orthopedic appliances, is contraindicated for use in patients with active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection. The system is also contraindicated for use in patients with known or suspected metal allergies.

Potential Adverse Events:
Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

Note: As with any surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, and permanent pain and/or deformity. Rarely, some complications may be fatal.

Warnings and Precautions:
The surgeon should be aware of the following when using metallic implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. The size and shape of human bones present limiting restrictions on the size and strength of implants. No implant can be expected to withstand the unsupported stresses of full weight bearing.
2. Correct handling of the implant is extremely important. Contouring the metal devices is to be avoided where possible. If contouring is necessary, it should not be bent sharply, reverse bent, notched or scratched. All of these operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.
3. Single use only. Reuse of devices labeled as single-use (e.g. implants, drills, tacks, trial rods) could result in injury or reoperation due to breakage or infection.
4. Non-sterile; the surgical mesh and instruments are sold non-sterile, and therefore, must be sterilized before each use.
5. Postoperative care is important. The patient should be instructed in the limitations of the metallic implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
6. All implants are intended for SINGLE USE ONLY. Any used implant should be discarded. Even though the device may appear undamaged, it may have small defects and internal stress patterns that may lead to fatigue failure.

MRI Compatibility Information:
The NGage Surgical Mesh System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of NGage Surgical Mesh System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Cleaning and Disinfection:
NGage Surgical Mesh System instruments and implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluid it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned after each use. Please refer to AW-10125289, Non-Sterile Surgical Instruments and Implants Reprocessing Guide, for instruction on the cleaning and disinfection of surgical instruments.

None of the instruments in the system require disassembly prior to cleaning.

Sterilization:
The NGage Surgical Mesh System components are supplied NON-STERILE. Please refer to AW-10125289, Non-Sterile Surgical Instruments and Implants Reprocessing Guide, for instruction on the sterilization of implants, components, and instrumentation.

Packaging:
Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The NGage Surgical Mesh System instruments and implants are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

Product Complaints:
Any Healthcare Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix US LLC, 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.

Further Information:
A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:
The implants, instruments and/or packaging material for the NGage Surgical Mesh System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

R_X Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician		
	See Instructions for Use		Manufacturer
	Orthofix.com/IFU	LOT	Lot Number
	Single Use Only Do Not Reuse	REF	Catalogue Number
	Provided Non-Sterile		

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