


93991 Galaxy Sterilization Box, Empty



LOWER EXTREMITY TRAY

Can accommodate:


Item	Part #	Description	Qty	Item	Part #	Description	Qty
1	93010	Large Clamp	14	7	932200	Rod D 12mm L 200mm	1
2	93020	Multi Screw Clamp	2	8	932150	Rod D 12mm L 150mm	1
3	932400	Rod D 12mm L 400mm	1	9	932100	Rod D 12mm L 100mm	1
4	932350	Rod D 12mm L 350mm	1	10	932030	Semi Circular Rod D 12mm large	1
5	932300	Rod D 12mm L 300mm	1	11	932020	Semi Circular Rod d 12mm medium	1
6	932250	Rod D 12mm L 250mm	1	12	932010	Semi Circular Rod d 12mm small	1

ORTHOFIX Srl
Via Delle Nazioni 9
37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000
Fax +39 045 6719380



Pins (all AO quick connect, all self drilling) and Clamp Interaction
Blue Clamp (12mm Rod) Pins, 6mm Shank only; 6mm & 4mm threads

93991 Galaxy Sterilization Box, Empty



UPPER EXTREMITY TRAY

Can accommodate:

Item	Part #	Description	Qty	Item	Part #	Description	Qty
1	93110	Medium Clamp	7	8	936200	Rod D 6mm L 200mm	1
2	93310	Small Clamp	7	9	936180	Rod D 6mm L 180mm	1
3	939300	Rod D 9mm L 300mm	1	10	936160	Rod D 6mm L 160mm	1
4	939250	Rod D 9mm L 250mm	1	11	936140	Rod D 6mm L 140mm	1
5	939200	Rod D 9mm L 200mm	1	12	936120	Rod D 6mm L 120mm	1
6	939150	Rod D 9mm L 150mm	1	13	936100	Rod D 6mm L 100mm	1
7	939100	Rod D 9mm L 100mm	1	14	936080	Rod D 6mm L 80mm	1
				15	936060	Rod D 6mm L 60mm	1

ORTHOFIX Srl
Via Delle Nazioni 9
37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000
Fax +39 045 6719380




Pins (all AO quick connect, all self drilling) and Clamp Interaction

Yellow Clamp (9mm Rods) Pins, 6mm Shank only; 6mm & 4mm threads

Green Clamp (6mm Rods) Pins, 4mm Shank only; 3mm threads

93992 Galaxy Instruments Sterilization Box, Empty



INSTRUMENTS FOR LOWER AND UPPER TRAYS

Can accommodate:

Item	Part #	Description	Qty	Item	Part #	Description	Qty
1	19940	Multi Screw Clamp Guide	1	10	19990	Drill Guide d 2.7mm	1
2	11138	Drill Guide D 4.8mm	2	11	1-1355001	Drill Bit D 2.7x127mm Coated - Quick Connect	1
3	11137	Screw Guide 80mm	2	12	19965	Tapered Trocar	1
4	1-1100201	Drill Bit D 4.8x240mm Coated - Quick Connect	2	13	93170	3mm Thread Screws Shaft Connection	1
5	11106	Drill Guide d 3.2mm	2	14	93150	Racheting T Handle	1
6	11102	Screw Guide 60mm	2	15	93155	Screw Shaft Connection	1
7	1-1300301	Drill Bit D 3.2x140mm Coated - Quick Connect	2	16	30017	Allen Wrench 5mm	1
8	19955	Trocar	1	17	93017	Wrench 5mm Shaft Connection	1
9	19985	Wrist Guide Template with Handle	1	18	93165	Screw Shaft Connector Quick Connect	1

Large Pin Tray

94990 6mm THREAD XCALIBER CYLINDRICAL SCREW TRAY, EMPTY

Can accommodate:


UPPER TRAY

Part #	Description	Qty
942630	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 260/30 QC	6
942640	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 260/40 QC	6
942650	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 260/50 QC	8
942660	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 260/60 QC	6
942670	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 260/70 QC	6


LOWER TRAY

Part #	Description	Qty
1-93050	Transfix Pin 50mm QC	2
1-93080	Transfix Pin 80mm QC	2
941630	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 180/30 QC	6
941640	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 180/40 QC	6
941650	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 180/50 QC	8
941660	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 180/60 QC	6
941670	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 6mm L 180/70 QC	6

UPPER TRAY



LOWER TRAY



Small Pin Tray

94991 4mm and 3mm THREAD XCALIBER CYLINDRICAL SCREW TRAY, EMPTY

Can accommodate:


UPPER TRAY

Part #	Description	Qty
947320	Self Drilling Xcaliber Cylindrical Screw Shaft D.4mm Thread 3mm L 100/20 QC	4
947325	Self Drilling Xcaliber Cylindrical Screw Shaft D.4mm Thread 3mm L 100/25 QC	4
948320	Self Drilling Xcaliber Cylindrical Screw Shaft D.4mm Thread 3mm L 120/20 QC	4
948325	Self Drilling Xcaliber Cylindrical Screw Shaft D.4mm Thread 3mm L 120/25 QC	4
948335	Self Drilling Xcaliber Cylindrical Screw Shaft D.4mm Thread 3mm L 120/35 QC	4
946420	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 180/20 QC	6
946430	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 180/30 QC	6
946440	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 180/40 QC	6


LOWER TRAY

Part #	Description	Qty
943420	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 100/20 QC	6
943430	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 100/30 QC	6
943440	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 100/40 QC	6
944420	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 120/20 QC	6
944430	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 120/30 QC	6
944440	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 120/40 QC	6
945420	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 150/20 QC	6
945430	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 150/30 QC	6
945440	Self Drilling Xcaliber Cylindrical Screw Shaft D.6mm Thread 4mm L 150/40 QC	6


UPPER TRAY



LOWER TRAY



93997 Galaxy Elbow Sterilization Box, Empty




ELBOW LOWER TRAY

Can accommodate:

Item	Part #	Description	Qty	Item	Part #	Description	Qty
1	93010	Large Clamp	1	10	11138	Drill Guide D 4.8mm	2
2	93020	Multi Screw Clamp	1	11	11137	Screw Guide 80mm	1
3	93110	Medium Clamp	2	12	11116	Drill Guide D 3.2mm L 80mm	2
4	93410	Elbow Hinge	1	13	19950	Drill Guide D 3.2mm L 100mm	2
5	932200	Rod D 12mm L 200mm	1	14	11102	Screw Guide 60mm	1
6	939150	Rod D 9mm L 150mm	1	15	1-1300301	Drill Bit D 3.2x140mm Coated - Quick Connect	2
7	30017	Allen Wrench 5mm	1	16	1-1100301	Drill Bit D 3.2x200mm Coated - Quick Connect	2
8	19940	Multi Screw Clamp Guide	1	17	11146	X-Wire without olive Ø 2 mm Length 150 mm	2
9	1-1100201	Drill Bit D 4.8x240mm Coated - Quick Connect	2	18	19955	Trocar	1

ORTHOFIX Srl
Via Delle Nazioni 9
37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000
Fax +39 045 6719380



**Pins (all AO quick connect, all self drilling)
and Clamp Interaction**

Blue Clamp (12mm Rod)

Pins, 6mm Shank only; 6mm & 4mm threads

Yellow Clamp (9mm Rods)


Pins, 6mm Shank only; 6mm & 4mm threads

ORTHOFIX[®]
ORTHOPEDICS

www.orthofix.com

GF-1306-PL-E0 D 01/22

93997 Galaxy Elbow Sterilization Box, Empty



ELBOW UPPER TRAY

Can accommodate:

Item	Part #	Description	Qty	Item	Part #	Description	Qty
1	932200	Rod D 12mm L 200mm	2	7	10025	Torque Wrench 6mm	1
2	93150	Racheting T Handle	1	8	93431	Humeral Distractor Clamp	1
3	93165	Screw Shaft Quick Connect	1	9	93432	Ulnar Distractor Clamp	1
4	93017	Wrench 5mm Shaft Connection	1	10	93420	Elbow Motion Unit	1
5	10017	Allen Wrench 6mm	1	11	93155	Screw Shaft Connection	1
6	93440	Wrench 5mm	1				

93998 Galaxy Shoulder Sterilization Box, Empty

SHOULDER IMPLANTS & INSTRUMENTS

Can accommodate:

Item	Part #	Description	Qty
1	93310	Small Clamp	4
2	93620	Wire Locking Clamp	3
3	936080	Rod D 6mm L 80mm	1
3	936100	Rod D 6mm L 100mm	1
3	936120	Rod D 6mm L 120mm	1
3	936140	Rod D 6mm L 140mm	1
4	30017	Allen Wrench 5mm	1
5	19975	Wire Targeting Device	1
6	19970	Wire Guide	2
7	19980	Wire Bender	1
8	91150	Bone Screw T Wrench	1
9	81031	Open End Wrench	1

Note:

99-93100 Wires – 2 per sterile pack
Not in tray, order separately.



Pins (all AO quick connect, all self drilling) and Clamp Interaction
Green Clamp (6mm Rods)
Pins, 4mm Shank only; 3mm threads

