


Minimize your incision,
Maximize your **visibility**

JANUS™
MIDLINE FIXATION SCREW




FIREBIRD® SPINAL FIXATION SYSTEM
PHOENIX® MINIMALLY INVASIVE SPINAL FIXATION SYSTEM

JANUS™

MIDLINE FIXATION SCREW

The JANUS Midline Fixation Screw is designed to achieve more cortical bone purchase in the medial to lateral trajectory when compared to traditional pedicle screws. With its modular screw design, maximized visualization is possible when working in a smaller midline incision or minimally invasive surgical approach.



Increased cortical bone purchase
43% increase in cortical region thread density*

**Compatible with Firebird® Spinal Fixation System
and Phoenix® Minimally Invasive Spinal Fixation System**


Self-drilling, self-tapping screws allow for fewer surgical steps

Threads extend to the tip for controlled insertion and tactile feedback


JANUS Midline Fixation Screw is available in both non-cannulated and cannulated versions.

*compared to Firebird 5.5mm x 30mm length pedicle screw

Non-Cannulated



Cannulated



IMPLANTS

Non-Cannulated Implant Caddy

Part #	Description
36-2425	4.5mm x 25mm
36-2430	4.5mm x 30mm
36-2435	4.5mm x 35mm
36-2440	4.5mm x 40mm
36-2445	4.5mm x 45mm
36-2525	5.5mm x 25mm
36-2530	5.5mm x 30mm
36-2535	5.5mm x 35mm
36-2540	5.5mm x 40mm
36-2545	5.5mm x 45mm
36-2625	6.5mm x 25mm
36-2630	6.5mm x 30mm
36-2635	6.5mm x 35mm
36-2640	6.5mm x 40mm
36-2645	6.5mm x 45mm

Optional Non-Cannulated Implants

Part #	Description
36-2450	4.5mm x 50mm
36-2455	4.5mm x 55mm
36-2460	4.5mm x 60mm
36-2465	4.5mm x 65mm
36-2550	5.5mm x 50mm
36-2555	5.5mm x 55mm
36-2560	5.5mm x 60mm
36-2565	5.5mm x 65mm
36-2650	6.5mm x 50mm
36-2655	6.5mm x 55mm
36-2660	6.5mm x 60mm
36-2665	6.5mm x 65mm
36-2725	7.5mm x 25mm
36-2730	7.5mm x 30mm
36-2735	7.5mm x 35mm
36-2740	7.5mm x 40mm
36-2745	7.5mm x 45mm
36-2750	7.5mm x 50mm
36-2755	7.5mm x 55mm
36-2760	7.5mm x 60mm
36-2765	7.5mm x 65mm

Cannulated Implant Caddy





Part #	Description
36-3425	4.5mm x 25mm
36-3430	4.5mm x 30mm
36-3435	4.5mm x 35mm
36-3440	4.5mm x 40mm
36-3445	4.5mm x 45mm
36-3525	5.5mm x 25mm
36-3530	5.5mm x 30mm
36-3535	5.5mm x 35mm
36-3540	5.5mm x 40mm
36-3545	5.5mm x 45mm
36-3625	6.5mm x 25mm
36-3630	6.5mm x 30mm
36-3635	6.5mm x 35mm
36-3640	6.5mm x 40mm
36-3645	6.5mm x 45mm

Optional Cannulated Implants





Part #	Description
36-3450	4.5mm x 50mm
36-3455	4.5mm x 55mm
36-3460	4.5mm x 60mm
36-3465	4.5mm x 65mm
36-3550	5.5mm x 50mm
36-3555	5.5mm x 55mm
36-3560	5.5mm x 60mm
36-3565	5.5mm x 65mm
36-3650	6.5mm x 50mm
36-3655	6.5mm x 55mm
36-3660	6.5mm x 60mm
36-3665	6.5mm x 65mm
36-3725	7.5mm x 25mm
36-3730	7.5mm x 30mm
36-3735	7.5mm x 35mm
36-3740	7.5mm x 40mm
36-3745	7.5mm x 45mm
36-3750	7.5mm x 50mm
36-3755	7.5mm x 55mm
36-3760	7.5mm x 60mm
36-3765	7.5mm x 65mm

INSTRUMENTS





Drills

	Part #	Description
	36-1124	4.5mm
	36-1125	5.5mm
	36-1126	6.5mm
	36-1127	7.5mm


Non-cannulated Cortical Bone Tap, Modular

	Part #	Description
	36-1224	4.5mm
	36-1225	5.5mm
	36-1226	6.5mm
	36-1227	7.5mm

Cannulated Cortical Bone Tap, Modular



	Part #	Description
	36-1424	4.5mm
	36-1425	5.5mm
	36-1426	6.5mm
	36-1427	7.5mm

Handle

	Part #	Description
	69-1030	Modular Handle

ADDITIONAL INSTRUMENTS

Drill Guide

	Part #	Description
	36-1012	Drill Guide
	36-1834	Midline Screw Driver

36-7200 Instruments



Top Level

Bottom Level

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience.

 Orthofix
3451 Plano Parkway
Lewisville, Texas 75056-9453 USA
1.214.937.3199
1.888.298.5700
www.orthofix.com

 Medical Device Safety Services (MDSS):
Schiffgraben 41
30175, Hannover
Germany
+49 511 6262 8630
www.mdss.com

Rx Only
CE²⁷⁹⁷

Orthofix products or services referenced herein are trademarks or registered trademarks of Orthofix Medical Inc. and its group of companies. Any rights not expressly granted herein are reserved.