



LUMBOSACRAL FIXATION


Unity[®]

Lumbosacral Fixation System






Robust keyed-features prevent implant rotation during placement, minimizing potential for great vessel trauma



Offset on plate holder offers unobstructed access to bone screw holes

Unity LX

Smooth surface and low profile design help reduce potential vascular trauma



Easy-to-read directional arrows help ensure proper implant alignment



Optional anterior or anterolateral placement for intraoperative flexibility



Lateral-to-medial screw orientation helps reduce risk of foraminal encroachment



The Unity LX™ Lumbar Fixation Plate



MATCHING THE RIGHT DEVICE TO THE ANATOMY EASILY, WITH CERTAINTY.

IMPLANTS—UNITY 51

43-1017*	Unity 51 Lumbosacral Fixation Base Plate	17mm
43-1019	Unity 51 Lumbosacral Fixation Base Plate	19mm
43-1021	Unity 51 Lumbosacral Fixation Base Plate	21mm
43-1023	Unity 51 Lumbosacral Fixation Base Plate	23mm
43-1025	Unity 51 Lumbosacral Fixation Base Plate	25mm
43-1028	Unity 51 Lumbosacral Fixation Base Plate	28mm
43-2017*	Unity 51 Lumbosacral Fixation Cover Plate	17mm
43-2019	Unity 51 Lumbosacral Fixation Cover Plate	19mm
43-2021	Unity 51 Lumbosacral Fixation Cover Plate	21mm
43-2023	Unity 51 Lumbosacral Fixation Cover Plate	23mm
43-2025	Unity 51 Lumbosacral Fixation Cover Plate	25mm
43-2028	Unity 51 Lumbosacral Fixation Cover Plate	28mm
43-6525	Bone Screw, Primary	25mm
43-6530	Bone Screw, Primary	30mm
43-6535	Bone Screw, Primary	35mm
43-7525	Bone Screw, Rescue	25mm
43-7530	Bone Screw, Rescue	30mm
43-7535	Bone Screw, Rescue	35mm
43-0104	Tack	
43-1090	Unity System Case (Empty)	

IMPLANTS—UNITY LX

43-3019	Unity LX Lumbar Fixation Base Plate	19mm
43-3021	Unity LX Lumbar Fixation Base Plate	21mm
43-3023	Unity LX Lumbar Fixation Base Plate	23mm
43-3025	Unity LX Lumbar Fixation Base Plate	25mm
43-3028	Unity LX Lumbar Fixation Base Plate	28mm
43-3031	Unity LX Lumbar Fixation Base Plate	31mm
43-4019	Unity LX Lumbar Fixation Cover Plate	19mm
43-4021	Unity LX Lumbar Fixation Cover Plate	21mm
43-4023	Unity LX Lumbar Fixation Cover Plate	23mm
43-4025	Unity LX Lumbar Fixation Cover Plate	25mm
43-4028	Unity LX Lumbar Fixation Cover Plate	28mm
43-4031	Unity LX Lumbar Fixation Cover Plate	31mm
43-6525	Bone Screw, Primary	25mm
43-6530	Bone Screw, Primary	30mm
43-6535	Bone Screw, Primary	35mm
43-7525	Bone Screw, Rescue	25mm
43-7530	Bone Screw, Rescue	30mm
43-7535	Bone Screw, Rescue	35mm
43-0104	Tack	
43-1091	Unity LX Lumbar Fixation Plate Case (Empty)	

INSTRUMENTATION

43-0099	Base Plate Holder, Long
43-0100	Base Plate Holder, Short
43-0103	Tack Holder / Driver
43-0104	Tack
43-0105	Fixed Guide
43-0106	Drill
43-0107	Tap
43-0108	Probe

43-0109	Awl
43-0111	Cover Plate Holder / Driver, Self-Retaining
43-0112	Hex Driver, Self-Retaining, Long
43-0113	Hex Driver, Self-Retaining, Short
43-1065	Gray Torque Limiting T-Handle, 50 in-lbs
55-1030	Ratcheting Handle, Straight
55-1031*	Ratcheting Handle, T-Style
74-0050	Mallet

* available upon request