


REEF™ TO

POSTERIOR INTERBODY SYSTEM
SALES BROCHURE

REEF TOPOGRAPHY™


Undercut macrostructures designed to promote bony interlocking^{5,+}



UP TO 22% MORE

NanoMetalene® Endplate Surface Area*

Endplate features place graft material in direct contact with the endplates.



UP TO 41% MORE

NanoMetalene Aperture Surface Area*



Aperture features secure graft within the aperture during interbody placement.

*Data on file.

NANOMETALENE® TECHNOLOGY

Submicron titanium layer molecularly bonded to entire PEEK implant



Titanium surfacing resists wear debris.¹



Radiolucent for post-op fusion assessment.²

Mechanical properties of PEEK unaltered, providing stiffness on par with bone.¹

Preclinical results show greater bone ongrowth on NanoMetalene® vs. PEEK^{3,†}



Rough submicron topography encourages integration.^{3,4}

REEF™ TO

The SeaSpine® Reef™ TO (TLIF Oblique) Interbody System is designed for posterior lumbar interbody procedures and accommodates both insert and rotate and direct impact insertion techniques. Implants feature NanoMetalene® surface technology and Reef Topography™ for greater titanium surface area and improved biomechanical stability while maintaining the physical properties of PEEK.^{1,3,†}

Interbody Features

- Reef Topography and NanoMetalene surface technology
- Footprint options: 9 x 24mm, 9 x 28mm, 11 x 28mm and 11 x 32mm
- Height options: 7–16mm*
- Lordotic options: 5°, 10°, 15°**



Instrumentation to allow for both open and minimally invasive techniques



Offset handle



T-handle

*Note: Min./Max. heights vary by footprint.


**TLIF implants offer 30° offset lordosis for oblique placement

PRECLINICAL EVALUATION



Preclinical Evaluation of Reef Topography™

Clinically relevant endplate-sparing sheep interbody fusion model results comparing NanoMetalene® (NM) implants with and without Reef Topography.

Endplate Undercut Macrostructures




Aperture Undercut Macrostructures



Undercut Macrostructures results in

~3x increase in mechanical stability^{5,†}



REEF TOPOGRAPHY MAY DRIVE EARLIER AND IMPROVED BIOMECHANICAL STABILITY^{5,†}

FUSION ENGINEERED™

Deliberate Design. Driven by Science.

SeaSpine® is dedicated to providing a comprehensive portfolio of innovative, procedurally-focused products strategically designed to work together to drive fusion. The latest advancements in bone biology and materials science guide the development of in-house manufactured advanced orthobiologics and proprietary spinal implant technology engineered to address the many nuances of spinal pathology. SeaSpine products can be tailored to meet individual patient needs, delivering both clinical and economic value to patients, surgeons, and hospital systems.



COMPLEMENTARY TECHNOLOGIES

The combination of SeaSpine's spinal implant technologies and orthobiologics capitalizes on both the design and the science, unleashing a favorable environment for fusion.



PROCEDURAL VERSATILITY

SeaSpine's comprehensive and procedurally-focused product portfolio offers flexibility that allows a tailored approach to each unique surgical situation.



SCIENTIFICALLY DRIVEN

SeaSpine is committed to relying on data to drive each facet of our technologies toward the ultimate goal of fusion.

¹Preclinical testing, such as animal studies, may not be indicative of human results.

²Results from mechanical testing. Data on file. TR-0010-11-01


³Results from imaging study. Data on file. TR-0010-11-01

⁴Walsh, et al. The in vivo response to a novel Ti coating compared with polyether ether ketone: evaluation of the periphery and inner surfaces of an implant. Spine Journal 2018 Jul; 18(7): 1231-1240

⁵NanoMetalene SEM images on file. TR-0094-19-01

⁶Results from preclinical in vivo testing. Data on File. D0003269

For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com | seaspine.com

 **SeaSpine Orthopedics Corporation**
5770 Armada Drive
Carlsbad, CA 92008 USA
TEL 760.727.8399 USA | FAX 760.727.8809



Outside USA

TEL +1.760.727.8399 | FAX +1.760.727.8809

INTERNATIONAL INQUIRIES intlcustomer@seaspine.com

ALL OTHER INQUIRIES customerservice@seaspine.com

SeaSpine, the SeaSpine logo, Accell, NanoMetalene and OsteoStrand are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. Reef Topography, Reef and WaveForm are trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2021 SeaSpine Orthopedics Corporation. All rights reserved. D0003915B 01