



ADMIRAL™ ACP

ANTERIOR CERVICAL FIXATION SYSTEM
SALES BROCHURE

Admiral™ ACP is a comprehensive and complete anterior cervical plating system designed to strike the optimal balance between strength, profile and construct rigidity. Its specialized instrumentation combined with unique plate features make Admiral ACP an efficient and reproducible ACDF solution.




90° turn locks screws securely in place




Large graft window


Comprehensive Implant Offerings



- Available in 4.0mm and 4.3mm diameter, 10–18mm lengths
- Fixed and variable, self-tapping and self-drilling options



No added profile from screw angulation or locking mechanism



OPTIMIZED INSTRUMENTATION FOR IMPROVED WORKFLOW



Osteophyter

The Osteophyter's unique cutting surface matches the plate's curvature to aid in ventral bony preparation.



Trial Drill Guides

Trial Drill Guides allow for precise pilot hole placement to ensure shortest plate possible.



Threaded Driver

Novel Threaded Driver provides robust control during screw insertion.



Temporary Fixation Pins

Temporary Fixation Pins can be inserted through the screw hole or locking cover.

PROCEDURAL SOLUTIONS


DIFFERENTIATED AND COMPLEMENTARY TECHNOLOGIES

Admiral™ ACP

ANTERIOR CERVICAL
FIXATION SYSTEM

WaveForm® C

Anterior Cervical Interbody



OsteoStrand® Plus

100% Demineralized Bone Fibers with
Benefits of Accell® Bone Matrix



For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com | seaspine.com

 **SeaSpine Orthopedics Corporation**
5770 Armada Drive
Carlsbad, CA 92008 USA
TEL 760.727.8399 USA | FAX 760.727.8809



Outside USA

TEL +1.760.727.8399 | FAX +1.760.727.8809
INTERNATIONAL INQUIRIES intlcustomer@seaspine.com
ALL OTHER INQUIRIES customerservice@seaspine.com

Warning: Applicable laws restrict these products to sale by or on the order of a physician. SeaSpine, the SeaSpine logo, Accell, OsteoStrand and WaveForm are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. Admiral is a trademark of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2021 SeaSpine Orthopedics Corporation. All rights reserved. D0004444B 01