




ANTERIOR DISC PREP

DISC PREP INSTRUMENTATION

ANTDISCPREP2 TRAY CONFIGURATION

Top Level – Spreader Inserter and Tips, Pituitaries, Kerrisons, Bone Spreader and Bone Rongeur



Spreader Inserter and Modular Tips




Spreader Inserter Assembled




ANTDISCPREP2 TRAY CONFIGURATION

Middle Level – Cobbs, Curettes and Rasp




Bottom Level – Curettes




ANTDISCPREP2 TRAY CONFIGURATION

Rake Curette



Single Sided Rasp



ANTDISCPREP2 Catalog

Part Number	Part Description
95-6990	ALIF Spreader
95-6191	ALIF Modular Spreader Inserter
95-6192	ALIF Modular Distractor Tip – 8mm
95-6193	ALIF Modular Distractor Tip – 9mm
95-6194	ALIF Modular Distractor Tip – 10mm
95-6195	ALIF Modular Distractor Tip – 11mm
95-6196	ALIF Modular Distractor Tip – 12mm
95-6197	ALIF Modular Distractor Tip – 13mm
95-6970	Cobb, 18mm
95-6971	Cobb, 23mm
95-6972	Cobb, 28mm
06-40-0258	Anterior Bone Rongeur
95-6974	Straight Cup Curette, Medium
95-6975	Straight Cup Curette, Large
95-6976	Straight Cup Curette, Extra Large

Part Number	Part Description
95-6977	Up Biting Curette, Small
95-6979	Up Biting Curette, Large
95-6980	Up Biting Currette, Extra Large
95-6982	Down Biting Curette, Medium
95-6985	Rake Curette
95-6986	Single Sided Rasp
93-5556	Pituitary, 13 in., 5mm
93-5569	Kerrison Punch, 6mm
93-5566	Kerrison Punch, 4mm
93-5568	Kerrison Punch, 2mm
95-6989	Scalpel Handle #3, 12 in.
95-6987	Rake curette – Angled 15°
95-6990*	ALIF Spreader

*Optional Instrument

For more information or to place an order, please contact:
 TEL 866.942.8698 | FAX 877.558.6227
 customerservice@seaspine.com | seaspine.com

Outside USA
 TEL +1.760.727.8399 | FAX +1.760.727.8809
 INTERNATIONAL INQUIRIES intlcustomer@seaspine.com
 ALL OTHER INQUIRIES customerservice@seaspine.com

