



Daytona[®] Small Stature



Daytona[®] Small Stature

Pedicle Screw Fixation System


Daytona[®] Small Stature Pedicle Screw Fixation System

DESIGN RATIONALE

The Daytona[®] Small Stature Pedicle Screw Fixation System is designed to address standard to complex deformity cases in patients who need a lower profile construct due to smaller anatomy constraints. It combines uniquely designed implants and materials with innovative instrumentation to create a highly flexible and intuitive system for surgeons and patients.

Indicated for Adolescent Idiopathic Scoliosis (AIS)

SYSTEM FEATURES




Hooks



Crossbars



Axial Connectors




CoCr and Ti Locking Caps



Rod-to-Rod Connectors

CoCr TEXTURED ROD

- 4.5mm diameter
- Straight and precontoured





Lines for dorsal identification and thoracolumbar junction



Daytona Small Stature Screw vs Daytona Screw
~20% reduced profile design

Cannulated Polyaxial Screws



6.5mm Diameter


5.5mm Diameter

4.5mm Diameter

7.5mm Diameter

Double Lead Thread

ADDITIONAL COMPLIMENTARY SOLUTIONS



Daytona[®]
Pedicle Screw
Fixation System



Ventura[™]
NanoMetalene[®]
Posterior Lumbar Interbody



Cabo[™] ACP
Anterior Cervical Fixation System



For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com

Outside USA
TEL + 1.760.727.8399 | FAX + 1.760.727.8809
INTERNATIONAL INQUIRIES intlcustomer@seaspine.com
ALL OTHER INQUIRIES customerservice@seaspine.com

seaspine.com



SeaSpine Orthopedics Corporation
5770 Armada Drive
Carlsbad, CA 92008 USA
TEL 760.727.8399 USA | FAX 760.727.8809

Warning: Applicable laws restrict these products to sale by or on the order of a physician.

SeaSpine, the SeaSpine logo, Daytona, and NanoMetalene are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. Cabo, Malibu, and Ventura are trademarks of SeaSpine Orthopedics Corporation or its subsidiaries. ©2018 SeaSpine Orthopedics Corporation. All rights reserved. Printed in USA. D00001105D 01