



# ISOTIS™ CERVICAL FACET FUSION

FACET PREP INSTRUMENTS AND STRAND® CERVICAL PELLETS  
SALES BROCHURE

**IsoTis™ Cervical Facet Fusion** is a novel procedural solution for reproducible biologic placement in the facet joint and lateral mass screw pathway preparation. Facet Prep Instruments and Strand® Cervical Pellets are packaged together within one product box to deliver a complete fusion kit, intended for use in conjunction with posterior cervical fixation systems. Featuring a set of single-use instruments and precision-shaped demineralized bone fiber pellets, IsoTis Cervical Facet Fusion enables:



**ACCURATE & CONTROLLED**

lateral mass screw preparation



**MAXIMIZED**

decorticated surface area within the facet joint




**SIMPLIFIED**

biologic delivery with precision-shaped demineralized bone fiber pellets



**EFFICIENT & REPRODUCIBLE**


method for robust facet fusion



# ISOTIS™ CERVICAL FACET FUSION


FACET PREP INSTRUMENTS AND STRAND® CERVICAL PELLETS

## FACET PREP INSTRUMENTS



### Cannulated Facet Drill Guide


accesses the facet joint and positions cervical screw pilot holes



### 5.0mm Facet Drill

creates a precise defect for facet joint preparation

**Lateral Mass Drill**  
(2.25mm) for 3.5mm screws, with adjustable Depth Sleeves



### Bone Graft Tamp

delivers biologic into the facet joint



## STRAND® CERVICAL PELLETS

- 100% demineralized bone fibers for maximized osteoinductive potential
- Entangled fibers provide osteoconductive scaffold
- Gentle expansion upon hydration *in situ* enables complete facet defect fill for fusion
- Available in 4- or 10-packs




# ISOTIS™ CERVICAL FACET FUSION

FACET PREP INSTRUMENTS AND STRAND® CERVICAL PELLETS

## Compatibility with SeaSpine® posterior cervical fixation system

### NorthStar® OCT


Posterior Cervical Fixation



Polyaxial screws  
featuring 90°  
of angulation



Strand® Cervical Pellets replaced with metal  
dowels for demonstration purposes



For more information or to place an order, please contact:  
TEL 866.942.8698 | FAX 877.558.6227  
customerservice@seaspine.com | [seaspine.com](http://seaspine.com)

#### Outside USA

TEL +1.760.727.8399 | FAX +1.760.727.8809  
INTERNATIONAL INQUIRIES [intlcustomer@seaspine.com](mailto:intlcustomer@seaspine.com)  
ALL OTHER INQUIRIES [customerservice@seaspine.com](mailto:customerservice@seaspine.com)

#### ISOtis OrthoBiologics, Inc.

2 Goodyear, Irvine CA 92618  
TEL 800.550.7155 | FAX 800.471.3248 | WEB [seaspine.com](http://seaspine.com)  
IsoTis OrthoBiologics, Inc. is a member of the SeaSpine  
Orthopedics Corporation family of companies.  
Made in the USA

