




SEASpine[®] JAZZ[™] BANDS

POSTERIOR BAND SYSTEM
SALES BROCHURE

BEST-IN-CLASS INTEGRATION

The SeaSpine® JAZZ™ Bands are temporary orthopedic implants intended to provide short-term stabilization for levels T1-L5 as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures.



JAZZ is the most recognized name in sublaminar banding with the widest range of implants on the market. Integration with SeaSpine posterior fixation systems, such as Mariner®, Malibu™, and Daytona®, provide more procedural options that are both streamlined and agile for the surgeon.



The combination of a pedicle screw and a sublaminar band can be used for top of construct reinforcement, an adjunct to fusion, or for pediatric deformity correction. The application of both is shown to support a stiffer construct than a pedicle screw alone.

SEASPINE® JAZZ™ BANDS

POSTERIOR BAND SYSTEM

JAZZ™ BAND SYSTEM OPTIONS



JAZZ™ BAND

Malleable stainless steel strip helps to facilitate passage of the band around bony structures



JAZZ CONNECTOR

- Secures band to a Ti or CoCr rod by a single locking screw
- Torque applied to screw transfers compressive force to band by way of rod
- Smooth connector design protects fraying of band



JAZZ PASSER BAND

Passer band available with multiple passer instrument angles designed to help surgeons navigate around bony structures

JAZZ LOCK

- Rodless zip-tie fixation system unique to the market
- Quick and easy one step crimping mechanism
- Optimized instrumentation for a short learning curve

SEASPINE® JAZZ™ BANDS

POSTERIOR BAND SYSTEM



For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com | seaspine.com

 **SeaSpine Orthopedics Corporation**
5770 Armada Drive
Carlsbad, CA 92008 USA
TEL 760.727.8399 USA | FAX 760.727.8809

Outside USA
TEL +1.760.727.8399 | FAX +1.760.727.8809
INTERNATIONAL INQUIRIES intlcustomer@seaspine.com
ALL OTHER INQUIRIES customerservice@seaspine.com

Warning: Applicable laws restrict these products to sale by or on the order of a physician.
SeaSpine, the SeaSpine logo, Daytona, and Mariner are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. JAZZ and Malibu are trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2019 SeaSpine Orthopedics Corporation. All rights reserved. Printed in USA. D0002838B 01