




Mariner®




Mariner®
Strength in Versatility

Mariner[®] Posterior Fixation System

PRODUCT POSITION

The Mariner Posterior Fixation System is a premier pedicle screw system featuring modular threaded technology, with thoughtfully designed instrumentation, packaged in a streamlined offering. Mariner is designed to provide surgeons with **Strength in Versatility** to facilitate posterior lumbar fixation.




SYSTEM FEATURES

- Modular Threaded Technology
- Six Unique Head Options
- Motion-Limiting Technology
- 5.5 and 6.0mm Rod Versatility
- User-Friendly Rod Reduction Applications
- Concise Three Tray System
- Dual-Diameter Shank Technology



Deformity Head

- Providing motion-limiting strength medial or lateral for techniques
- Max angulation for optimal rod capture
- 65° in three planes



T30 Drive




Large Rhomboid Threads

- Minimizing cross-threading
- Maximizing implant performance

Trauma Heads

- Providing motion-limiting strength cephalad or caudal for superior distraction techniques
- 65° in three planes



REDUCTION OPTIONS



Rod Rocker

- Providing 10mm of reduction
- Cantilever Rod into Inferior or Superior Screw



Pistol Persuader

- Providing 30mm of reduction
- Controlled sequential reduction



Threaded Reducer

- Providing 30mm of reduction
- Apply additional force to seat rod



For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com

Outside USA
TEL + 1.760.727.8399 | FAX + 1.760.727.8809
INTERNATIONAL INQUIRIES intlcustomer@seaspine.com
ALL OTHER INQUIRIES customerservice@seaspine.com

seaspine.com



SeaSpine Orthopedics Corporation

5770 Armada Drive
Carlsbad, CA 92008 USA
TEL 760.727.8399 USA | FAX 760.727.8809

Warning: Applicable laws restrict these products to sale by or on the order of a physician.

SeaSpine, the SeaSpine logo and Mariner are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2018 SeaSpine Orthopedics Corporation. All rights reserved. Printed in USA. D0001175D 01