

POSTERIOR DISC PREP


DISC PREP & DECOMPRESSION INSTRUMENTATION

PDISCPREP TRAY CONFIGURATION

Top Level – Pituitaries & Handles





Middle Level – Currettes & Rasps



PDISCPREP TRAY CONFIGURATION

Bottom Level – Distractors & Scrapers



Rotary Scraper*

*Cutting edges on both sides



Paddle Distractor

Key Instrument Working Ends



Angled Curette, #3


Straight Curette, #3

Teardrop Curette, 8mm


Rasp, 8mm

PDECOMPRESS TRAY CONFIGURATION

Top Level – Curettes, Kerrisons, Chisels, Funnel



Bottom Level – Soft Tissue Retractors



For more information or to place an order, please contact:
TEL 866.942.8698 | FAX 877.558.6227
customerservice@seaspine.com | seaspine.com



Warning: Applicable laws restrict these products to sale by or on the order of a physician. SeaSpine and the SeaSpine logo are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. ©2021 SeaSpine Orthopedics Corporation. All rights reserved. D0001710B 01