


# MINIMALLY INVASIVE ACCESS SYSTEMS

SALES BROCHURE

## THE SEASPINE® TUBE RETRACTOR SYSTEM

The SeaSpine Tube Retractor System is designed as a percutaneous access port for minimally invasive posterior lumbar decompression and interbody procedures. Multiple working tube dimensions, bright surgical field illumination, and secure table stabilization complements the SeaSpine portfolio of minimally invasive posterior lumbar products.



# MINIMALLY INVASIVE ACCESS SYSTEMS

## IPASSAGE™ RETRACTOR


The iPassage™ Retractor adjusts to a variety of anatomical and surgical conditions. It offers a main retractor frame, a third blade option, and a counter retractor allowing for third and fourth blade attachment. Blade and hook options add flexibility and the ability to angulate up to 30° increases potential visibility to three spinal levels.

### Open Retractor Architecture

Optimizing surgical field exposure

### Secure Retractor Base

Beveled and radiolucent fan blades to minimize soft tissue obstruction



### Control Muscle Creep

Independent blade adjustment for precise adaptation to anatomy

### Optimal Visualization

Reduces shadows and dissipates heat through the retractor blades

# MINIMALLY INVASIVE ACCESS SYSTEMS


## TUBE RETRACTOR

Item	Size
Dilators	7, 13, 18, 22mm
Radiolucent Tube Diameter	18 and 22mm
Radiolucent Tube Length	40–120mm

## IPASSAGE™ RETRACTOR


Item	Size
Dilators	8, 13, 18, 22mm
Base Retractors (Semi-tubular Radiolucent Blades)	18 and 22mm (Closed Diameter) 50–120mm (Lengths)
Counter Retractor and Center Blade (Straight Hooks, Straight and Fan Blades)	40–120mm (Lengths)

For more information or to place an order, please contact:  
TEL 866.942.8698 | FAX 877.558.6227  
customerservice@seaspine.com | [seaspine.com](http://seaspine.com)

Tube Retractor:  
 **SeaSpine Orthopedics Corporation**  
5770 Armada Drive  
Carlsbad, CA 92008 USA  
TEL 760.727.8399 USA | FAX 760.727.8809



**Outside USA**  
TEL + 1.760.727.8399 | FAX + 1.760.727.8809  
INTERNATIONAL INQUIRIES intlcustomer@seaspine.com  
ALL OTHER INQUIRIES customerservice@seaspine.com

iPassage:  
 **Tedan Surgical Innovations LLC**  
12615 W Airport Blvd Suite 200  
Sugar Land, TX 77478 USA  
TEL 713.726.0886

Warning: Applicable laws restrict these products to sale by or on the order of a physician. SeaSpine and the SeaSpine logo are registered trademarks of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries. iPassage is a trademark of SeaSpine Orthopedics Corporation or its subsidiaries in the United States and/or other countries.  
©2020 SeaSpine Orthopedics Corporation. All rights reserved. Printed in USA. D0002982B 01