

Flex TLSO SP

DCT-64 The Flex TLSO SP immobilizes and restricts mid and lower spine movement (T9- S1). Ideal for post-operative recovery, the Flex TLSO SP limits forward movement during the healing period. A lightweight and highly supportive frame equipped with our tool-free TX2 ActivAdjust Sternal Kit, breathe through mesh belt, and moisture wicking lining increase wearability and patient compliance.

Superior Post-Operative Support

Limits forward movement during the healing period to optimize recovery



Highlights

- Provides immobilization and motion restriction
- Patent-pending dual pulley system telescopes the thoracic core independently from the lumbar brace; this revolutionary mechanism secures a neutral spine alignment, and an increased low load assisted thoracic stretch
- Removeable sternal pad
- Over/under design

Indications

- Conditions requiring gross immobilization of the thoracic and lumbar spine
- Post-surgical immobilization
- Kyphosis
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures

PDAC Approved: L0464