

Trauma Collar

DCT-0180 The Trauma Collar is a revolutionary multi-post cervical collar that provides a customized contoured fit and consistent reapplication while minimizing skin-related issues and patient error. Designed to address the most common collar issues, the Trauma Collar features an advanced height adjustment system that allows for 60 unique sizes, a quick-clip system that delivers consistent reapplication, and advanced open-cell padding that increases air flow and decreases pressure.

Patented

Precision Dial

Allows for a smooth, straightforward fitting experience improving compliance



Highlights

- Height setting from 2 7/8" to 6 3/8"
- 60 possible sizing options
- Customizable front and back panels
- Breathable open-cell padding
- Sizing landmarks (chin and straps)
- Moisture-wicking, latex-free

Indications

- Cervical spine immobilization
- Cervical stenosis
- Post-trauma wear
- Arthritis
- Cervical disc herniation
- Cervical radiculopathy
- Post-operative care

Consistent Reapplication

Quick-clip closure system and sizing landmarks eliminate patient error when reapplying the brace in supine or ambulatory patients

PDAC Approved: L0180