

Pro Wrist

DCT-3916 The Pro Wrist is a multipurpose upper extremity orthosis that provides hand and wrist support. The Pro Wrist is a progressive pain management solution for pre-, post-, and non-operative hand and wrist issues, including carpal tunnel syndrome.



Highlights

- Immobilizes and restricts hand/ wrist movement
- Provides support and stabilization
- Reduces median nerve pressure
- Can be set in flexion, extension, or full range of motion

Indications

- Carpal tunnel syndrome
- Tendonitis
- Wrist immobilization
- Distal radius fracture
- Mild to moderate sprains and strains
- Post cast healing
- Rheumatoid arthritis

Precision ROM Dial

Allows providers ultimate control in setting flexion or extension to optimize patient recovery and reduce median nerve pressure

Maximum Comfort

Quick-lace compression design and moisture wicking neoprene sleeve maximize patient compliance

PDAC Approved: L3915/L3916